

Database Merging

1 Backup database

1.1 Backup old database

- Make a safety copy of the database you want to import

1.2 Backup current database


- Utilities→Configuration→General→Storage→ Enter Password: EOPTM
- **A:** Write down path and name of current database
- →Press Backup...
- Store database to an external memory
- *Default database backup name: CopyOfDatabase*

2 Import a old database

- Utilities→Configuration→General→Storage→ Enter Password: EOPTM → Import...
- Select the database you want to import
- Press OK (database will be imported)
- Database merging is finished

If this error message occurs, proceed with step 3

Note: An old database format applies, it need to be converted



3 Converting old database

- Utilities→Configuration→General→Storage→ Enter Password: EOPTM → Select...
- Select database you want to open (the one you want to merge with the existing)
- Press OK (database will be imported and converted into the new format, EasyOne Pro does restart)
- Utilities→Configuration→General→Storage→ Enter Password: EOPTM → Import...
- Browse previous data base, path noted under **A:**
- Change name of backup database: Old database converted