

EasyOne FlowTube

Breathing Tube for Infection Control

Avoid cross-contamination and reduce cleaning to an absolute minimum with the single-patient-use EasyOne FlowTube.



All parts exposed to patient's breath are **single-use** and **disposable**

No expiration date

No special storage conditions

Technical data

Material

EasyOne FlowTube: Polyethylene (PP)
Bag: Polyethylene(PE)

Transport and storage information

Temp. -20 to 50°C
Rel. Humidity 5 to 90%
Athmosph. Pressure 500 to 1060 hPa
no expiration date

MDR (EU) 2017/745 // MDD 93/42/EEC

Class IIa

Dimension (single EasyOne FlowTube)

21 x 38 x 150 mm (H x W x D), 8 g

Compatible with

EasyOne Air

Order information

Part Number

Description

5050-50

EasyOne FlowTube, individually wrapped, dispenser box of 50 pcs

5050-200

EasyOne FlowTube, individually wrapped, box of 200 pcs (dispenser box included)