

Spirette™

Breathing Tube for Infection Control

Avoid cross-contamination and reduce cleaning to an absolute minimum with the single-patient-use Spirette.



All parts exposed to patient's breath are **single-use** and **disposable**

No expiration date

No special storage conditions

Technical data

Material

Spirette: Linear low density polyethylene (LLDPE)
Bag: Polyethylene (PE)

Transport and storage information

Temp. -20 to 50°C
Rel. Humidity 5 to 95%
Atmosph. Pressure 500 to 1060 hPa
no expiration date

MDR (EU) 2017/745 // MDD 93/42/EEC

Class IIa

Dimension (single Spirette)

25 x 38 x 140 mm (H x W x D), 8 g

Compatible with

EasyOne, Easy on-PC, EasyOne Pro and EasyOne Pro LAB

Order information

Part Number

Description

2050-1

Spirette, individually wrapped, dispenser box of 50 pcs

2050-5

Spirette, individually wrapped, box of 200 pcs (dispenser box included)

2050-10

Spirette, individually wrapped, box of 500 pcs (without dispenser box)