

EasyOne FlowTube

Breathing mouthpiece for Infection Control



n d d **EasyOne FlowTube** mouthpieces avoid cross-contamination between patients and reduce cleaning of the device to an absolute minimum.

The **EasyOne FlowTube** mouthpieces have been engineered to fit perfectly into n d d's device to guarantee reliable and high quality measurement results.

All parts exposed to patient's breath are **single-use and disposable**

No expiration date

Technical Data

Material	EasyOne FlowTube: Polypropylene (PP) Bag: Polyethylene (PE)
Transport and storage information	Temperature: -25 to 70°C Rel. Humidity: 5 to 90%, non-condensing Atmospheric Pressure: 400 to 1060 hPa
Dimension (single EasyOne FlowTube)	21 x 38 x 150 mm (H x W x D), 8 g
Compatible with	EasyOne Air

Order Information

Part number:	Product:
5050-50	EasyOne FlowTube, individually wrapped, dispenser box of 50 pcs
5050-200	EasyOne FlowTube, individually wrapped, box of 200 pcs (dispenser box included)
5050-500	EasyOne FlowTube, individually wrapped, box of 500 pcs (without dispenser box)