

# EasyOne FlowTube

Breathing Tube for  
Infection Control

**Avoid cross-contamination  
and reduce cleaning to an  
absolute minimum with  
the single-patient-use  
EasyOne FlowTube.**



All parts exposed to patient's breath  
are **single-use** and **disposable**

No expiration date

No special storage conditions

## Technical data

### Material

EasyOne FlowTube: Polyethylene (PP)  
Bag: Polyethylene (PE)

### Transport and storage information

Temp. -20 to 50°C  
Rel. Humidity 5 to 90%  
Athmosph. Pressure 500 to 1060 hPa  
no expiration date

### MDR (EU) 2017/745 // MDD 93/42/EEC

Class IIa

### Dimension (single EasyOne FlowTube)

21 x 38 x 150 mm (H x W x D), 8 g

### Compatible with

EasyOne Air

## Order information

### Part Number

### Description

5050-50	EasyOne FlowTube, individually wrapped, dispenser box of 50 pcs
5050-200	EasyOne FlowTube, individually wrapped, box of 200 pcs (dispenser box included)
5050-500	EasyOne FlowTube, individually wrapped, box of 500 pcs (without dispenser box)